



V1.0.6

KINOVA SDK USER GUIDE

DISCLAIMER

Kinova, Kinova SDK, and Kinova's logo are trademarks of Kinova Inc., herein referred to as Kinova. All other brand and product names are trademarks or registered trademarks of their respective corporations.

The mention of any product does not constitute an endorsement by Kinova. This manual is furnished under a lease agreement and may only be copied or used within accordance with the terms of such lease agreement. Except as permitted by such lease agreement, no part of this publication may be reproduced, stored in any retrieval system, or transmitted, in any form or by any means, electronic, mechanical, recording, or otherwise, without prior written consent of Kinova.

The content of this manual is furnished for informational use only, is subject to change without notice, and should not be construed as a commitment by Kinova. Kinova assumes no responsibility or liability for any errors or inaccuracies that may appear in this document.

Changes are periodically made to the information herein; these changes will be incorporated into new editions of this publication. Kinova may make improvements and/or changes in the products and/or software programs described in this publication at any time.

Address any questions or comments concerning this document, the information it contains or the product it describes to:

- support@kinovarobotics.com

Kinova may use or distribute whatever information you supply in any way it believes appropriate without incurring any obligations to you.

Copyright © 2013-2015 Kinova. All rights reserved.

TABLE OF CONTENTS

DISCLAIMER..... I

TABLE OF CONTENTS II

OVERVIEW 1

Content.....1

INSTALLATION..... 2

Windows 72

Windows 85

Ubuntu 12.04/14.04 installation7

DEVELOPMENT CENTER..... 10

General settings10

Advanced settings14

Monitoring.....18

Virtual joystick21

Trajectory planner25

Resources.....27

Examples.....28

CONTACTING SUPPORT 30

OVERVIEW

The Kinova SDK is a complete set of interface, documentation, examples and software tools that help the developer interact with any Kinova product. It is available under Ubuntu and Windows systems.

Content

- The development center
- A set of project examples
- User guide
- HTML documentation on all APIs
- Tools to configure your product


INSTALLATION

Windows 7


NEW INSTALLATION

If you have any Kinova products already installed on your computer, please refer to the [UPDATE](#) section. If it is your first installation of a Kinova product, follow the procedure below.

1. Download and install [Microsoft Visual C++ Redistributable\(x86, x64\)](#)
2. Execute the installer named **KinovaSDKInstaller**.
3. Connect the robot to your computer via USB.
4. Power on the robot.
5. In the Windows Control Panel, open the device manager and wait for the computer to detect a Custom USB Device.



6. Right click on the Kinova product that appeared and install the drivers that were copied on your disk when you executed the KinovaSDKInstaller. Assuming that you've installed the SDK in the default folder, it should be located at C:\Program Files (x86)\KinovaSDK\
7. A Windows Security window may appear:




8. Choose to install the driver software anyway.
9. The SDK is ready!

UPDATE

This section explains how to install the Kinova SDK if you have any Kinova products already installed on your computer. This is necessary because the Kinova SDK has a new driver that is different from the one used in previous versions of Kinova software products like Jacosoft. Note that it is possible to have different driver versions on different USB ports. As an example, you could install the previous USB driver on USB port A and when the robot is connected on that port, Jacosoft is available. At the same time, you could install the new driver on USB port B and when your robot is connected to that port, the new SDK is available.

The procedure below explains how to uninstall a driver completely before installing the new Kinova SDK.


1. Connect the robot to your computer via USB.
2. Power on the robot.
3. In Windows Control Panel, open the device manager and look for the Custom USB Device called Jaco Robotical Arm.



4. Right click on it and choose uninstall



5. Check the option Delete the driver software for this device and click OK.
6. Stay in the device manager and right click on your computer. Choose Scan for hardware changes.




7. From here, follow the procedure in the section [NEW INSTALLATION](#).


Windows 8.1

NEW INSTALLATION


1. Move your mouse to the bottom-right corner of the screen and click on the Settings button.




2. Click on the Change PC settings button located at the bottom of the menu.



3. Select the sub-menu Update and recovery.
4. Select the sub-menu recovery.
5. Click on the Restart Now button located in the Advanced startup section.
6. Wait until it restarts.
7. Select Troubleshoot.



8. Select Advanced options.
9. Select Startup Settings.
10. Click on the Restart button.
11. Press the F7 key on your keyboard to disable driver signature enforcement.



12. Right click on the bottom-left corner of the screen and select Device Manager.
13. From there, follow the same procedure as for [Windows 7](#).


UPDATE

This section explains how to install the Kinova SDK if you have any Kinova products already installed on your computer. This is necessary because the Kinova SDK has a new driver that is different from the one used in previous versions of Kinova software products like Jacosoft. Note that it is possible to have different driver versions on different USB ports. As an example, you could install the previous USB driver on USB port A and when the robot is connected on that port, Jacosoft is available. At the same time, you could install the new driver on USB port B and when your robot is connected to that port, the new SDK is available.


To install the Kinova SDK, follow the same [driver uninstallation](#) procedure as in Windows 7 and proceed afterward to Windows 8.1 [new installation](#).

Ubuntu 12.04/14.04 installation


1. Execute the bash script named *installSDK32.sh* if your computer's architecture is 32 bits or *installSDK64.sh* if your computer's architecture is 64 bits.
2. The script will ask for root permission to install a package named kinova-api.
3. The Kinova SDK Installer Setup window will appear to install the remainder of the SDK.



4. Choose an installation folder (default: /opt/kinova)



5. Choose which components you need to install.




6. Click on the next button until installation is completed.
7. The SDK is installed. To launch the Development Center, execute [installation folder]/GUI/KinovaSDK.sh

DEVELOPMENT CENTER

General settings

SETTINGS



① Advanced HOME / RETRACTED position

- Activate

Check to enable the advanced HOME / RETRACTED feature.
- List of positions

Contains the list of positions (angular or Cartesian) that represents the new path that the robot will follow between the HOME and the RETRACTED position.
- Add

Click to add the current robot's position at the end of the

positions list.

Remove	Click to remove the selected position from the position list.
Insert below	Click to add the current robot's position below the selected position in the positions list.
Insert above	Click to add the current robot's position above the selected position in the positions list.
Clear all	Click to clear the entire positions list.

② Options


Retracted angle	The angle of the default RETRACTED position shown by the image aside.
Handedness	Select the radio button to set the robot in right handed mode or in left handed mode.
Invert fingers	Check to invert the closing finger when a 3 fingers robot is closing only 2 fingers. This mode is mainly used with a left handed configuration.
Control mode	Displays the current control mode. It can be either angular or Cartesian.
Switch button	Select to toggle between the angular control mode and the Cartesian control mode. Note that if the robot's position cannot switch to Cartesian mode because of a singularity, the robot will stay in angular mode.
Velocity	Select the Cartesian velocity of the robot.

- Activate offset

Check to apply a position offset on the end effector. The offset is described by the text fields X, Y and Z.
- Save

Press to send the modification to the robot. If you perform a modification in the Option panel nothing will be applied unless the save button is pressed. The exception is the velocity which is applied directly when modified.

UPDATE



1 Firmware update

- White panel

Displays the complete path to the HEX file that will be uploaded to the robot when the button update is pressed.

Browse	Press to choose a HEX file on your disk that contains a firmware update.
Update	Press to update your robot with the chosen HEX file.
Restore factory settings	Restore the robot with its factory settings.

② Versions

Main firmware	Displays the main firmware's code version.
Actuator X	Displays the code's version of the actuator X where X is a number greater than 0.
Finger X	Displays the code's version of the finger X where X is a number greater than 0.

Advanced settings

ROBOT 1

Advanced settings

Advanced

PID

P: 0

Command filter: 0

Address: 16

Apply

I: 0

Derivative filter: 0

Apply to all

D: 0

Error filter: 0

Zero position

Address: 0

Apply

Apply to all

Address

current: 0

new: 0

Apply

Torque zero

Address: 0

Apply

Apply to all

Info

Please note that the default address of the actuator #1 which is the actuator near the base is set to 16 and you increment by one until you have reach the last finger. Example: actuator #4's address is 19.

Reference frame

☐ Fixed frame

☒ Rotating frame

Reactive force control

☒ Activated


☐ Not activated

Current limitation

☒ Activated

☐ NotActivated

1 PID (Proportional, Integral and derivative)



P

The proportional part of the controller. (Suggested range [0, 2])

I	The integral part of the controller.
D	The derivative part of the Controller. (Suggested range [0, 0.1])
Command filter	Filters the commands sent to the actuator. The filter is of first order and the parameter value is a frequency in rad/s. The suggested range is [0, 500].
Derivative filter	Filters the error derivative signal (derivative of (desired position - measured position)). The filter is of first order and the parameter value is a frequency in rad/s. The suggested range is [0, 500].
Error filter	Filters the error signal (desired position - measured position). The filter is of first order and the parameter value is a frequency in rad/s. The suggested range is [0, 500].



Setting these parameters out of the suggested range can severely damage the robot and it is not covered by the warranty.

② Zero position

Address	The address of the actuator to modify.
Apply	Apply the <i>zero position</i> . The actual position of the targeted actuator will now be 180°.
Apply to all	Apply the <i>zero position</i> to every actuator of the robot.



Setting these parameters incorrectly can severely damage the robot and it is not covered by the warranty.

③ Address

Current The current address of the actuator.

New The new address to assign.

Apply Apply the new address on the actuator.



Setting these parameters incorrectly can severely damage the robot and it is not covered by the warranty.

④ Torque zero

Address The address of the actuator on which the *torque zero* will be applied if you press Apply.

Apply Apply the *torque zero*. The actual torque of the targeted actuator will now be considered as 0

Apply to all Apply the zero position on every actuator of the robot.



Setting these parameters incorrectly can severely damage the robot and it is not covered by the warranty.

⑤ Reference frame

Fixed frame Set the robot to fixed frame. When a translation is performed, the orientation stays the same.

Rotating frame Set to rotating. When a translation is performed, the orientation follows the first actuator.

⑥ Reactive force control

Activated Select to activate the reactive force control.

Not activated Select to deactivate the reactive force control.

⑦ Current limitation

Activated Select it to activate the current limitation.


Not activated Select it to deactivate the current limitation.



Setting these parameters incorrectly can severely damage the robot and it is not covered by the warranty.

Monitoring

GENERAL



Monitoring ROBOT 1

General

System

Running time 20467 seconds

- ▶ Power & Sensors
- ▶ Status
- ▶ Limitation
- ▶ Control

Peripheral list

- ▶ Comm interface
- ▶ Rotary actuator
- ▶ Rotary actuator
- ▶ Rotary actuator
- ▶ Rotary actuator
- ▶ Rotary actuator
- ▶ Rotary actuator
- ▶ Linear actuator
- ▶ Linear actuator
- ▶ Linear actuator
- ▶ Joystick

KINOVAROBOTICS.COM

① **System**

Displays general information stored in the robot.

② **Peripheral list**

Lists all peripherals detected on the robot's communication bus. Each peripheral can be expanded to display its type, address, port and code version.

ANGULAR

ROBOT 1

Monitoring

General

Angular

Cartesian

Angular

	Current	Torque	Command	Position	Encoder	Temp
	(A)	(N*m)	(°)	(°)	(°)	(°C)
Actuator 1	0.000	0.138	270.01	270.17	000.00	30.450
2	0.000	1.806	150.00	150.05	000.00	26.050
3	0.000	5.235	27.01	26.80	000.00	27.190
4	-0.006	0.062	268.00	267.89	000.00	28.930
5	0.000	0.040	5.00	5.18	000.00	31.170
6	0.000	-0.181	100.01	99.89	000.00	32.770
Finger 1	(A)		(u)	(u)	(u)	(°C)
	0.00	N/A	5412.00	5412.00	000.00	35.55
2	0.00	N/A	5412.00	5412.00	000.00	34.21
3	0.00	N/A	5406.00	5406.00	000.00	34.49

KINOVAROBOTICS.COM

① Angular

Displays all the information on every actuator. It includes the current, torque, command, position, encoder value, temperature, velocity, acceleration along the X axis, acceleration along the Y axis and acceleration along the Z axis.

CARTESIAN

ROBOT 1

Monitoring						
General		Cartesian				
		Current	Force	Command	Position	Velocity
		(A)	(N)	(m)	(m)	(m/s)
Angular	X	N/A	-0.998	0.0277	0.0279	0.00000
	Y	N/A	-1.392	-0.1700	-0.1689	0.00000
	Z	N/A	0.721	0.2991	0.2976	0.00000
Cartesian	Theta X	N/A	(N*m) 0.372	(Rad) 2.1456	(Rad) 2.1498	(Rad/s) 0.00000
	Theta Y	N/A	0.224	0.0755	0.0759	0.00000
	Theta Z	N/A	0.046	-0.1834	-0.1793	0.00000
	Finger 1	0.000	N/A	(u) 5412.00	(u) 5412.00	(u/s) 0000.00
	Finger 2	0.000	N/A	5412.00	5412.00	0000.00
	Finger 3	0.000	N/A	5406.00	5406.00	0000.00


KINOVAROBOTICS.COM

① Cartesian

Displays all the information on the end effector. It includes the force, command, position and velocity.

Virtual joystick

ANGULAR



① Angular control

Actuator X + Press and hold to move the actuator X counter clockwise where X is a number great than 0.

Actuator X - Press and hold to move the actuator X clockwise where X is a number great than 0.

② Gripper

Open 3 fingers Press and hold to open a 3 finger gripper.

Close 3 fingers Press and hold to close a 3 finger gripper.

Open 2 fingers Press and hold to open a 2 finger gripper.


Close 2 fingers Press and hold to close a 2 finger gripper.

③ Angular velocity

Text field Set the velocity of a specific actuator.

Save Send the new velocity configuration to the robot.

CARTESIAN



① Translation

- | | |
|-----------------|--|
| Button + | Press and hold to move the end effector along the positive axis described by the button. |
| Button - | Press and hold to move the end effector along the negative axis described by the button. |

② Orientation

- | | |
|-----------------|--|
| Button + | Press and hold to move the end effector around the axis described by the button counter clockwise. |
| Button - | Press and hold to move the end effector around the axis described by the button clockwise. |

③ Gripper


- | | |
|------------------------|---|
| Open 3 fingers | Press and hold to open a 3 finger gripper. |
| Close 3 fingers | Press and hold to close a 3 finger gripper. |
| Open 2 fingers | Press and hold to open a 2 finger gripper. |
| Close 2 fingers | Press and hold to close a 2 finger gripper. |

④ Cartesian velocity

- | | |
|-----------------------------|---|
| Translation velocity | Set the translation velocity of the end effector. |
|-----------------------------|---|

Orientation velocity Set the orientation velocity of the end effector.

KEYBOARD



① Control (Angular or Cartesian)

Combo box Select a key that will be mapped with the specified movement. Note that this panel can be changed by pressing the angular or Cartesian button.


② Gripper

Combo box Select a key that will be mapped with the specified gripper movement.

③ Others

- Move HOME** Select a key that will be mapped with the Move HOME functionality.
- Activate keyboard** Check to activate the keyboard control with the robot. Note that the robot's joystick will always have a higher control priority than the keyboard.
- Cartesian** Press to switch the control panel to Cartesian.
- Angular** Press to switch the control panel to Angular.

Trajectory planner



Trajectory ROBOT 1

Position

Actuator 1 deg

Actuator 2 deg deg/s

Actuator 3 deg

Actuator 4 deg

Actuator 5 deg deg/s

Actuator 6 deg

Finger 1

Finger 2

Finger 3

Options

Control mode: Angular

☐ Enable limitations ☐ Loop trajectory

☐ Enable fingers

Delay (sec)

Trajectory

KINOVAROBOTICS.COM

① Position (Angular or Cartesian)

Left section	Displays the current position or the position of a selected trajectory point.
Right section	Displays the max velocity of the robot for the selected trajectory point.

② Trajectory


Import	Press to import a trajectory. The trajectory must have been saved in a KTJ file on disk.
Export	Press to export a trajectory on disk in a KTJ file.
Position list	Displays the list of positions contained in the trajectory. Each point can be selected and their data will be displayed on the position panel.
Get	Press to display the current robot's position on the position panel.
Add	Press to add a point in the list based on the information on the position panel.
Quick add	Press to execute a combination of a Get and Add.
Update	Press to update the selected point in the list with the information from the position panel.
Reach	Press to move the robot to the selected point in the list.
Delete	Press to delete the selected point in the list.

Clear all	Press to clear all the positions in the list.
Stop	Press stop the robot movement. It will also cancel any trajectory loops.
Launch	Press execute the trajectory. If the Loop Trajectory check box is checked, the robot will execute the trajectory until the stop button is pressed.

③ Options

Move HOME	Select a key that will be mapped with the Move HOME functionality.
Activate keyboard	Check to activate the keyboard control with the robot. Note that the robot's joystick will always have a higher control priority than the keyboard.
Cartesian	Press to switch the control panel to Cartesian.
Angular	Press to switch the control panel to Angular.

Resources



These are shortcuts to support elements found on Kinova's website.

Examples

SOURCE

Angular control

Source

```
#include <iostream>
#include <dlfcn.h>
#include <vector>
#include "Lib_Examples/Kinova.API.CommLayerUbuntu.h"
#include "Lib_Examples/Kinova.API.UsbCommandLayerUbuntu.h"
#include "Lib_Examples/KinovaTypes.h"
#include <stdio.h>
#include <unistd.h>

using namespace std;

int main()
{
    int result;
    QuickStatus data;
    AngularPosition currentCommand;

    cout << "Angular control" << endl;

    //Handle for the library's command layer.
    void * commandLayer_handle;

    //Function pointers to the functions we need
    int (*MyInitAPI)();
    int (*MyCloseAPI)();
    int (*MySendBasicTrajectory)(TrajectoryPoint command);
    int (*MyGetDevices)(KinovaDevice devices[MAX_DEVICES], int &result);
    int (*MySetActiveDevice)(KinovaDevice device);
    int (*MyMoveHome)();
    int (*MyInitFingers)();
```

①


KINOVAROBOTICS.COM

① Source

Code

Displays the code that reflects the example.

CONSOLE



1 Console

Output Displays the output of the example's execution in a console like window.

CONTACTING SUPPORT

If you need help or have any questions about this product, this guide, or the information detailed in it, please contact a Kinova representative at:

- Support@KinovaRobotics.com

We value your comments!

To help us assist you more effectively with problem reports, the following information is required when contacting Kinova or your distributor support

- Product serial number
- Date/Time of the problem
- Environment where the problem occurred (per example 30° Celsius, raining, ...)
- Actions performed immediately before the problem occurred



6110 Doris-Lussier
Boisbriand, Quebec (Canada)
J7H 0E8

www.KinovaRobotics.com